

Tested disinfectants

The manufacturers of the following cleaning and disinfecting agents have carried out material suitability tests and confirmed compatibility with our brush strips. In some cases, the tests were carried out in high doses to simulate compatibility over a long period of time. All disinfectants should be dosed and used according to the manufacturer's recommendations.

Long-term use – permanent disinfection

Manufacturer	Tested disinfectants	Tested concentration in %	Duration of the test	Approval by manufacturer
	BTS 2000	1,8	90 days	
	BTS 2001			
	Des AM blau			
	Virodex			
	Calgonit DS 683	5	19 days	
	SaniPerfect VS61	2	70 days	
	Suredis VT1			
	P3 topax 990	3	28 days	
	Deptal AS	2	33 days	
		20		
	1+1 Wofasteril SC super	1	30 days	
	MENNO® Clean	6 and 8	120 days	
	MENNO® Florades			
	A-QUASAN®			
	Lerasept AC	2	30 days	
	Lerasept SOC	6		

Tested disinfectants

Short-term use – cleaning the brush strips

Manufacturer	Tested disinfectants	Tested concentration in %	Duration of the test	Approval by manufacturer
	Blaue Sandy	100	4 days	
	Goldene Sandy			
	ASiRAL SRC	100	10 days	
	Suma San con. D10.1 conc	3	30 days	
	P3 topax 66	5	28 days	
	Lerasept AC	2	30 days	
	Lerasept SOC	2,5		